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
One Sample, One GEM Well, One Flowcell


In this example you have one **sample** that is processed through one **GEM well** (a set of partitioned cells from a single 10x Chromium™ Chip channel) and sequenced on one **flowcell**. In this case you would generate FASTQs using `cellranger mkfastq`, and run `cellranger count` as described in [Single-Sample Analysis](#).

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One Sample, One GEM well, Multiple Flowcells


In this example you have one sample that is processed through one GEM well then you generate one library which is sequenced across multiple flowcells. This may be done to increase sequencing depth, for example. In this case all of the reads can be combined in a single instance of the `cellranger count` pipeline. This process is described in [Specifying Input Fastqs](#).

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One Sample, Multiple GEM Wells, One Flowcell

In this example you have one sample that is processed through multiple GEM wells. This is often done when conducting technical replicate experiments. The libraries from the GEM wells are then pooled onto one flowcell and sequenced. In this case you demultiplex the data from the sequencing run and then run the libraries from each GEM well through a separate instance of `cellranger count`. Once those are completed, you can perform a combined analysis using `cellranger aggr`, as described in [Multi-Library Aggregation](#). (See figure above.)

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